

LIST OF PLATES

Plate no.	Title	Page no.
4.3.3.1	Amber color glass bottles used for collection of fried oil samples	4.12
4.3.4.1.2.4	Preparation of column chromatograph for determination of total polar components in fried oils	4.15
4.3.4.1.2.6	Thin layer chromatograph sheet for separation of polar and non-polar components in fried oils	4.17
4.3.4.1.3.3	Determination of peroxide value by titration method	4.19
4.3.4.1.4.3.1	Dilution of oil sample with isooctane in 25 ml volumetric flask	4.21
4.3.4.1.4.3.2	Shimadzu UV-1201 used for determination of p-anisidine value	4.22
4.3.4.1.4.3.3	Addition of p-anisidine reagent in diluted samples of oils for second absorbance (As)	4.22
4.3.4.2.1	Abbe refractometer used for determination of refractive index of oil	4.29
4.3.4.2.2	Lovibond tintometer for estimation of color of oil used for frying french fries and bhajias	4.30
5.2.1.1	Panelists performing sensory evaluation	5.33
5.2.1.2	Puffiness/Absorption of oil in french fries at the end of 21 h of intermittent frying	5.33
5.2.2.1	Appearance of bhajias fried in CSO at initial hours of frying	5.41
5.2.2.2	Change in appearance of bhajias fried in CSO at the end of frying hours	5.41
5.3.1.1	Thin layer chromatograph shows the polar and non polar components separated by column chromatography from french fried CSO	5.69
5.3.1.2	French fries fried GNO and CSO fatty acid profile chromatograms obtained at 0 and 25 h of intermittent frying	5.72
5.3.1.4	GNO color before 25 h of intermittent french fries frying	5.74
5.3.1.5	GNO color after 25 h of intermittent french fries frying	5.74
5.2.3.1	Bhajias fried GNO and CSO fatty acid profile chromatograms obtained at 0 and 25 h of intermittent frying	5.84
5.3.2.2	CSO color before and after 25 h of intermittent bhajias frying	5.86
5.3.2.3	GNO color before and after 25 h of intermittent bhajias frying	5.86